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1. Bjarnason NH. & Christiansen, C. Early Response in Biochemical Markers Predicts Long-term Response in Bone Mass During Hormone Replacement Therapy in Early Postmenopausal Women. *Elsevier* 26; 561-569 (2000).
2. Bonde M, Garnero P, Fledelius C, Qvist P, Delmas PD, Christiansen C. Measurement of bone degradation products in serum using antibodies reactive with an isomerized form of an 8 amino acid sequence of the C-telopeptide of Type I collagen. *J Bone Miner Res* 12:1028-1034 (1997).
3. Bonde M, Qvist P, Fledelius C, Riis BJ, Christiansen C. Applications of an enzyme immunoassay for a new marker of bone resorption (CrossLaps®) - follow up on hormone replacement therapy and osteoporosis risk assessment. *J Clin Endocrinol&Metab* 80:864-8 (1995).
4. Buclin T, Rochat MC, Burckhardt P, Azria M, Attinger M. Bioavailability and biological efficacy of a new oral formulation of salmon calcitonin in healthy volunteers. *J Bone Miner Res* 17;8-1478 (2002).
5. Burgeson RE. New collagens, new concepts. *Annu Rev Cell Biol* 4:551-77 (1988).
6. Chailurkit LO, Aunphongpuwanart S, Ongphiphadhanakul B, Jongjaroenprasert W, Sae-Tung S, and Rajatanavin R. Efficacy of intermittent low dose alendronate in Thai Postmenopausal osteoporosis. *Endocr Res* 30:29-36 (2004).
7. Christgau S, Bitsch-Jensen O, Bjarnason NH, Gamwell Henriksen E, Qvist P, Alexandersen P, Bang Henriksen D. Serum CrossLaps for Monitoring the Response in Individuals Undergoing Antiresorptive Therapy. *Elsevier* 26; 505-511 (2000).
8. Christgau S, Rosenquist C, Alexandersen P, Bjarnason NH, Ravn P, Fledelius C, Herling C, Qvist P, Christiansen C. Clinical evaluation of the Serum CrossLaps® ELISA, a new assay measuring the serum concentration of bone-derive degradation products of type I collagen Ctelopeptides. *Clin Chem* 44:11;2290-2300 (1998).
9. Christiansen C, Tanko LB, Warming L, Moelgaard A, Christgau S, Qvist P, Baumann M, Wiczorek L, and Hoyle N. Dose dependent effects on bone resorption and formation of intermittently administered intravenous ibandronate. *Osteoporos Int* (2003).
10. Delmas PD. Markers of bone turnover for monitoring treatment of osteoporosis with antiresorptive drugs. *Osteoporos Int* 11;S66-S76 (2000).
11. Delmas PD, Hardy P, Dain MP. Monitoring Individual Response to Hormone Replacement Therapy with Bone Markers. *Elsevier* 26;553-560 (2000).
12. Fink E, Cormier C, Steinmetz P, Kindermans C, Bouc Le Y, Souberbielle JC. Differences in the Capacity of Several Biochemical Bone Markers to Assess High Bone Turnover in Early Menopause and Response to Alendronate Therapy. *Osteoporos Int* 11;295-303 (2000).
13. Garnero P, Gineys E, Riou JP, Delmas PD. Assessment of bone resorption with a new marker of collagen degradation in patients with metabolic bone disease. *J Clin Endocrinol&Metab* 79:780-5(19 94).
14. Pedersen BJ, Ravn P, Bonde M. Type I Collagen C-telopeptide Degradation products as Bone Resorption Markers. *J Clin Ligand Assay* 21;2:118-127 (1998).
15. Qvist P, Munk M, Hoyle N, Christiansen C. Serum and plasma fragments of C-telopeptides of type I collagen (CTX) are stable during storage at low temperatures for 3 years. *Clin Chim Acta* 350:167-173 (2004).
16. Ravn P, Neugebauer G, Christiansen C. Association between pharmacokinetics of oral ibandronate and clinical response in bone mass and bone turnover in women with postmenopausal osteoporosis. *Bone* 30;320-324 (2002).
17. Ravn P. & Clemmesen B. Biochemical Markers Can Predict the Response in Bone Mass During Alendronate Treatment in Early Postmenopausal Women. *Bone* 24, 237-244 (1997).
18. Reginster JY, Wilson KM, Dumont E, Bonvoisin B, Barrett J. Monthly Oral Ibandronate is Well Tolerated and Efficacious in Postmenopausal Women: Results From the Monthly Oral Pilot Study. *J Clin Endocrinol Metab* (2005).

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