COMP[®] ELISA

Enzyme Immunoassay

- Directions for Use, 3-12**
- Gebrauchsanweisung, 13-22**
- Mode d'emploi, 23-33**
- Instrucciones de uso, 35-44**
- Istruzioni per l'uso, 45-54**
- Aanwijzingen voor gebruik, 55-64**
- Instruções de Utilização, 65-74**
- Bruksanvisning, 75-84**
- Brugsanvisning, 85-94**

INTENDED USE

The COMP ELISA is a quantitative enzyme-linked immunosorbent assay for the determination of Cartilage Oligomeric Matrix Protein (COMP) in human serum. The presence of elevated levels, when considered in conjunction with other laboratory and clinical findings, is an aid in identifying aggressive destruction of joint tissue in diseases such as rheumatoid arthritis (RA).

SUMMARY AND EXPLANATION OF THE TEST

COMP was first described in 1992 by the research group of Professor Dick Heinegård (1), University of Lund, Sweden. The molecule is a pentameric protein of five identical disulfide-linked subunits (2) with a total molecular mass of 434 kDa (3). COMP was first described in cartilage (1) but has recently also been found in tendons (4) and synovial membranes (5). It has not been possible to detect COMP in culture medium from other connective tissue such as skin or lung tissue (6).

It is well documented (7,8,9,10) that when articular cartilage matrix is degraded by a disease process, protein fragments are produced and diffuse out into the joint fluid. Some of these proteins, such as COMP, subsequently appear in the bloodstream and can be used to monitor the progress of cartilage degradation in inflammatory joint diseases such as rheumatoid arthritis (8,11) and osteoarthritis (12,13, 14). A quantitative relation has been shown between the concentrations of COMP in serum and cartilage degradation (10) using radiographic changes as a surrogate clinical endpoint. This has been further proved in experimentally induced arthritis in animal models where the COMP serum level was highly correlated to the severity of arthritis (15,16) and also to the clinical joint score and histopathological signs of cartilage erosion (16).

PRINCIPLE OF THE PROCEDURE

COMP ELISA is a solid-phase, two-site enzyme immunoassay. It is based on the direct sandwich technique in which two monoclonal antibodies are directed against separate antigenic determinants on the COMP molecule. During incubation, COMP in the sample reacts with peroxidase-conjugated anti-COMP antibodies and anti-COMP antibodies bound to the microtitration well. A simple washing step removes unbound enzyme-labeled antibody. The bound conjugate is detected by a reaction with 3,3',5,5'-tetramethylbenzidine (TMB). This reaction is stopped by adding acid to give a colorimetric endpoint that is read spectrophotometrically.

REAGENTS

Each COMP ELISA kit contains reagents for 96 wells, which is sufficient for one calibrator curve, a blank, two controls and 40 samples in duplicate. For larger series of assays, use pooled reagents from packages bearing identical lot numbers. The expiry date for the complete kit is stated on the outer label. The recommended storage temperature is 2-8 °C.

Anti-COMP Coated Plate (Mouse monoclonal antibodies)	1 plate 96 wells 8-well strips	Ready for use Store at 2-8 °C until expiry date. For unused microtitration strips, reseal the bag and store at 2-8 °C for two months.
Calibrators (human COMP in buffer) Concentration indicated on vial label	5 vials 1.0 ml	Lyophilized Store at 2-8 °C until expiry date Add 1.0 ml distilled water per vial.* Reconstituted calibrator: store at 2-8 °C for 4 weeks or at -20 °C until expiry date.
Sample Buffer	2 vials 7.5 ml	Ready for use Store at 2-8 °C until expiry date.
Enzyme Conjugate 11x Peroxidase-anti-COMP (Mouse monoclonal antibodies, ~ 10 µg/ml)	1 vial 1.3 ml	Concentrate Store at 2-8 °C until expiry date. Preparation, see table. Do not freeze!
Conjugate Buffer Color coded blue	1 vial 13 ml	Ready for use Store at 2-8 °C until expiry date.
Wash tablets	2 tablets	Dissolve each tablet in 500 ml distilled water. Store wash solution at 2-8 °C. Use within 4 weeks.
Enzyme Substrate (TMB)	1 vial 22 ml	Ready for use Store at 2-8 °C until expiry date. Light sensitive!
Stop Solution** 0.5 M H ₂ SO ₄	1 vial 7 ml	Ready for use Store at 2-8 °C until expiry date.
Control 1 (human COMP in buffer)	1 vial 0.5 ml	Lyophilized Add 0.5 ml distilled water.* Reconstituted control: store at 2-8 °C for 4 weeks or at -20 °C until expiry date.
Control 2 (human COMP in buffer)	1 vial 0.5 ml	Lyophilized Add 0.5 ml distilled water.* Reconstituted control: store at 2-8 °C for 4 weeks or at -20 °C until expiry date.

*Note! Make sure that all lyophilized reagents are properly dissolved.

**Material Safety Data Sheet is available upon request

PREPARATION AND HANDLING OF SAMPLES AND ENZYME CONJUGATE

Samples

Blood should be collected by venipuncture, allowed to clot, and the serum separated by centrifugation. Specimens can be stored at room temperature (for shipping purpose only), at 2-8 °C if assayed within one week after collection, or at -20 °C if assayed later.

Heparin plasma can also be used.

Note! EDTA plasma and citrate plasma cannot be used.

Contact AnaMar for advice regarding measurement of COMP in other body fluids.

Preparation of samples and controls

Samples and controls should be diluted 1/10 in sample buffer.

(20 µl serum + 180 µl sample buffer.)

Preparation of Enzyme Conjugate

Prepare the enzyme conjugate by dilution of enzyme conjugate 11X, (1+10) in conjugate buffer.

Note: Conjugate volumes are preferably adjusted by diluting all at once. In other words, pour all of the buffer into the conjugate vial.

Mix gently. Store diluted conjugate at 2-8 °C for up to four weeks.

Or prepare the needed volume according to the table below. Mix gently.

Number of strips	Enzyme Conjugate 11X	Conjugate Buffer
4 strips	350 µl	3.5 ml
6 strips	500 µl	5.0 ml
12 strips	1 vial	1 vial

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use. Not for internal or external use in humans or animals.
- The contents of this kit and their residues must not be allowed to come into contact with ruminating animals or swine.

This kit contains reagents that may be infectious!

This kit contains reagents manufactured from human blood components. The source material has been tested by immunoassay for hepatitis B surface antigen and for antibodies for HIV virus and found to be negative. Nevertheless, all recommended precautions for handling blood derivatives should be observed. Please refer to HHS Publication no. (CDC) 88-8395 or corresponding local/national guidelines on laboratory safety procedures.

The Stop Solution consists of dilute sulfuric acid solution. Avoid exposure to bases, metals, or other compounds that may react with acids. Sulfuric acid is a poison and corrosive, which may be toxic if ingested. To prevent chemical burns avoid contact with skin and eyes.

PROCEDURE

Parameters of the Procedure

Volumes per well:

Sample	25 µl
Conjugate	100 µl
Enzyme Substrate	200 µl
Stop Solution	50 µl

Incubation Time and Temperature:

1st incubation	120 min. on a plate shaker at room temperature, 20-28 °C, (RT).
2nd incubation	15 min. at RT

Materials Required but not Provided

- 25 µl micropipette with disposable tips
- 50 µl, 100 µl and 200 µl repeating pipettes
- 1000 ml beaker
- Distilled water
- EIA plate reader with 450 nm filter
- Plate shaker
- Wash device for microtitration plates

Reference Material

COMP Calibrators are calibrated against an in-house standard correlated to the COMP inhibition assay described by Saxne and Heinegård in *Cartilage oligomeric matrix protein: A novel marker of cartilage turnover detectable in synovial fluid and blood*, published in British Journal of Rheumatology vol 31 pp 583-591 (1992).

Interference

Components have been added to reagents in the kit to prevent interference from mouse- and RF-antibodies present in the patient sera.

Test Procedure

Prepare a calibrator curve for each assay run. All reagents and samples must be brought to room temperature before use. Controls should be treated as samples. We recommend that each determination for calibrators and unknowns is performed in duplicate.

Add to anti-COMP wells

1. Calibrators	25 µl
2. Diluted controls and diluted samples	25 µl
3. Enzyme conjugate	100 µl
4. Incubate on a shaker and shake intensely (eg 250 rpm, Ø 16mm) for 2 hours at room temperature at setting.	
5. Wash 6 times with automatic washer, or: Aspirate the reaction volume. Add 350 µl Washing solution to each well. Aspirate completely. Repeat 5 times. After the final wash, invert and tap the plate firmly against absorbent paper.	
6. Add 200 µl enzyme substrate (TMB)	
7. Incubate for 15 minutes	
8. Add 50 µl stop solution.	
9. Put the plate on the shaker for approximately 5 seconds to ensure mixing of substrate and stop solution.	
10. Measure absorbance at 450 nm and evaluate.	

Quality Control

- The absorbance of the lowest COMP ELISA Calibrator should be greater than the absorbance of the COMP ELISA sample buffer.
- The COMP ELISA Control 1 and 2 concentrations should be within the range of the values stated on its label.

For each lot of controls, the concentrations of COMP have been determined using COMP ELISA from AnaMar. The mean value and standard deviation are stated on the respective vial label. Each mean value is based on six replicates in six consecutive runs.

As with all immunoassays, the results are affected by the testing procedures and equipment used by different laboratories. It is therefore recommended that each laboratory establish its own mean values and criteria of acceptance (e.g. 2 S.D. or 3 S.D.). A range of ± 2 S.D. should encompass 95% of the single determinations in a data set with normal distribution.


Example of Worksheet

Wells	Identity
1A-B	Sample buffer, (blank)
1C-D	Calibrators 0.4 U/L*
1E-F	Calibrators 0.7 U/L*
1G-H	Calibrators 1.2 U/L*
2A-B	Calibrators 1.8 U/L*
2C-D	Calibrators 3.2 U/L*
2E-F	Control 1
2G-H	Control 2
3A-B	Unknown1

*Concentration indicated on vial label.

Calibrator curve

A typical calibrator curve is shown below. Do not use this curve to determine actual assay results.

CALCULATION OF RESULTS

Computerized calculation

Computerized data reduction of absorbance for the Calibrators versus the concentration using Cubic Spline regression or four parameter regression may be performed to obtain the concentration of COMP. Multiply the concentration of the unknown samples with the dilution factor (eg. X 10).

Manual calculation


1. Plot the absorbance values obtained for the calibrators against the COMP concentration on a lin-log paper and construct a calibrator curve.
2. Read the concentration of the unknown samples from the calibrator curve.
3. Multiply the results by the dilution factor (x10).

EXPECTED VALUES

Expected values in blood donors

Good practice dictates that each laboratory establishes its own expected range of values. The following results may serve as a guide until the laboratory has gathered sufficient data of its own.

All samples were serum samples collected from blood donors in Sweden.

COMP in Early and Late RA

In a study of 183 rheumatoid arthritis patients, the concentration of COMP at inclusion was related to the risk of extensive small joint destruction after 5 years. Odds ratio per COMP unit in serum was 1.21 (confidence interval (CI) 1.07-1.36, $p < 0.002$). These data were used to estimate the relative risk of developing joint destruction above the 75th percentile (Larsen Score > 59 of 200) after 5 years (17). Other data published by Skoumal M *et al.* indicate that this relation to small joint destruction is also seen in patients with long-standing disease (18).

Interpretation

The following guidelines may be used to assess ongoing cartilage turnover and the risk of joint destruction in the future.

< 12 U/L	Lower risk of aggressive joint destruction*
12-15 U/L	Increasing risk of aggressive joint destruction*
> 15 U/L	High risk of aggressive joint destruction*

* An elevated ESR further increased the risk of aggressive joint destruction in RA (17).

PERFORMANCE CHARACTERISTICS

Detection limit⁽¹⁾

The detection limit is < 0.1 U/L.

Parallelism⁽¹⁾

Serum samples were diluted 10, 20, and 40 times using the sample diluent. The observed Obtained/Expected values for 1/20 and 1/40 were in the range 93-116%.

Hook Effect⁽¹⁾

Samples diluted 1/10 can have a concentration up to at least 200 U/L without giving incorrectly low concentrations.

Precision⁽¹⁾

Each sample was analysed in 4-replicates on six different occasions.

Sample	Mean value U/L	Coefficient of variation		
		within assay %	between assay %	total assay %
1	6.9	3.0	1.8	3.5
2	13.1	1.9	2.7	3.3
3	18.0	1.7	4.2	4.5

Notes

⁽¹⁾Studies performed at AnaMar Medical, Uppsala, Sweden

LIMITATIONS OF THE PROCEDURE

A definitive clinical diagnosis should not be based on any single diagnostic method, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

WARRANTY

The performance data presented here were obtained using the recommended procedure indicated. Any change or modification in the procedure not recommended by AnaMar Medical may affect the results, in which event AnaMar Medical disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and fitness for use.

In such an event, AnaMar Medical and its authorized distributors shall not be liable for damages indirect or consequential.

EXPLANATION OF SYMBOLS USED ON LABELS

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1. Hedbom E, Antonsson P, Hjerpe A, Aeschlimann, Paulsson M, Rosa-Pimentel E, *et al.* Cartilage matrix proteins. An acidic oligomeric protein (COMP) detected only in cartilage.
J Biol Chem 1992;267:6132-6.
2. Oldberg A, Antonsson P, Lindblom K, Heinegard D. COMP (cartilage oligomeric matrix protein) is structurally related to the thrombospondins.
J Biol Chem 1992;267:22346-50.
3. Zaia J, Boynton RE, McIntosh A, Marshak DR, Olsson H, Heinegard D *et al.* Post-translational modifications in cartilage oligomeric matrix protein. Characterization of the N-linked oligosaccharides by matrix-assisted laser desorption ionization time-of-flight mass spectrometry.
J Biol Chem 1997;272:14120-.
4. DiCesare P, Hauser N, Lehman D, Pasumarti S, Paulsson M. Cartilage oligomeric matrix protein (COMP) is an abundant component of tendon.
FEBS Lett 1994;354:237-40.
5. Di Cesare PE, Carlson CS, Stollerman ES, Chen FS, Leslie M, Perris R. Expression of cartilage oligomeric matrix protein by human synovium.
FEBS Lett 1997;412:249-52.
6. Recklies AD, Baillargeon L, White C. Regulation of cartilage oligomeric matrix protein synthesis in human synovial cells and articular chondrocytes.
Arthritis Rheum 1998;41:997-1006.
7. Neidhart M, Hauser N, Paulsson M, DiCesare PE, Michel BA, Hauselmann HJ. Small fragments of cartilage oligomeric matrix protein in synovial fluid and serum as markers for cartilage degradation.
Br J Rheumatol 1997;36:1151-60.
8. Månsson B, Carey D, Alini M, Ionescu M, Rosenberg LC, Poole AR, Heinegard D, Saxne T. Cartilage and bone metabolism in rheumatoid arthritis. Differences between rapid and slow progression of disease identified by serum markers of cartilage metabolism.
J Clin Invest 1995;95:1071-7.
9. Månsson B, Geborek P, Saxne T. Cartilage and bone macromolecules in knee joint synovial fluid in rheumatoid arthritis: relation to development of knee or hip joint destruction.
Ann Rheum Dis 1997;56:91-6.

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AnaMar Medical AB